

Item	Description
<b>Test Name</b>	Helix Comprehensive Cardiomyopathy Panel
<b>Test Type</b>	Cardio
<b>Catalog Number</b>	CCMP1
<b>Procedure Code</b>	H00223-2 (Helix)
<b>Test Description</b>	This panel evaluates 86 genes that have an established, primary association with cardiomyopathy, and several syndromic conditions of which cardiomyopathy is a feature.
<b>Genes Tested</b>	<i>ABCC9, ACAD9, ACADVL, ACTC1, ACTN2, AGL, ALMS1, ALPK3, BAG3, BMP10, BRAF, CDH2, CPT2, CRYAB, CSRP3, DES, DMD, DNAJC19, DOLK, DSC2, DSG2, DSP, DTNA, ELAC2, EMD, FHL1, FKRP, FKTN, FLNC, GAA, GLA, HCN4, HRAS, JPH2, JUP, KRAS, LAMP2, LMNA, LZTR1, MAP2K1, MAP2K2, MRAS, MTO1, MYBPC3, MYH7, MYL2, MYL3, MYLK3, MYPN, NEXN, NKX2-5, NRAS, PCCA, PCCB, PKP2, PLN, PPA2, PPCS, PRDM16, PRKAG2, PTPN11, RAF1, RBM20, RIT1, RYR2, SCN5A, SGCD, SHOC2, SLC22A5, SOS1, SOS2, SYNE2, TAFAZZIN, TBX20, TCAP, TMEM43, TMEM70, TNNC1, TNNI3, TNNI3K, TNNT2, TPM1, TRIM63, TTN, TTR, VCL</i>
<b>Genetics Information</b>	This test utilizes next-generation sequencing to detect single nucleotide variants, insertions and deletions up to 20 bp, and copy number variants in genes associated with hereditary forms of cardiomyopathy.
<b>Indications For Testing</b>	<p>Providing a genetic evaluation for individuals with a personal and/or family history suggestive of a hereditary form of cardiomyopathy.</p> <p>Establishing a diagnosis of a hereditary form of cardiomyopathy predisposition.</p>
<b>Clinical Descriptions</b>	<p>Cardiomyopathies are a broad spectrum of structural and functional disorders of the heart musculature. There are many different causes of cardiomyopathies, which range from environmental exposures to inherited genetic risk factors. In cases where an external cause is not identified, and/or a family history is suspicious for hereditary risk, diagnostic genetic testing may be indicated.</p> <p>Hypertrophic Cardiomyopathy (HCM) is a condition in which the heart muscle becomes abnormally thick. This leads to a decrease in the heart's ability to pump blood effectively. This can cause fatigue, shortness of breath, swelling of the legs, heart rhythm abnormalities and, in severe cases, heart failure or sudden cardiac arrest. Dilated Cardiomyopathy (DCM) is a condition in which the heart chambers become enlarged, weakening the heart muscle. In individuals with left ventricular noncompaction (LVNC), LVNC is characterized by endomyocardial trabeculations which can have variable effects on heart musculature, including ventricular dilation. Individuals with DCM or LVNC may be asymptomatic, or may be symptomatic with arrhythmia, left ventricular dysfunction, thromboembolic disease, and/or life-threatening arrhythmias. Arrhythmogenic cardiomyopathy (ACM) is characterized by fibrofatty infiltration of ventricle musculature which may present as arrhythmogenic right ventricular cardiomyopathy/dysplasia (ARVC) or arrhythmogenic left ventricular cardiomyopathy/dysplasia (ALVC). These findings are associated with increased risk for ventricular dysfunction, and arrhythmogenic events which may result in syncope or, in rare cases, cardiac arrest or sudden death.</p>

Item	Description
<b>Clinical Descriptions (continued)</b>	<p>It is important to note that in some cases HCM, DCM, LVNC and ARVC may be a feature of a larger syndromic condition. Hereditary forms of these conditions may follow autosomal dominant, autosomal recessive, X-linked or mitochondrial inheritance patterns. This panel does not assess mitochondrial inheritance. Note that some of these genes may also be associated with other unrelated conditions; this means that when undergoing this test, there is a possibility of incidentally detecting carrier status for, or predisposition to, one of these unrelated conditions.</p>
<b>Conditions</b>	<p>Hypertrophic cardiomyopathy (HCM)  Dilated cardiomyopathy (DCM)  Left ventricular non-compaction (LVNC)  Arrhythmogenic right ventricular cardiomyopathy (ARVC)  Transthyretin amyloidosis (ATTR-CM)  Brugada syndrome (BrS)</p> <p>In addition, one or more of the genes on this panel are associated with other conditions for which a predisposition to, or carrier status of, may incidentally be identified:</p> <p>Autosomal recessive Acyl-CoA dehydrogenase 9 deficiency (ACAD9)  Autosomal recessive Very long chain acyl-CoA dehydrogenase deficiency (ACADVL)  Autosomal recessive Glycogen storage disease type III (AGL)  Autosomal recessive Alstrom syndrome (ALMS1)  Autosomal dominant Cardiofaciocutaneous syndrome (BRAF, KRAS, MAP2K1, MAP2K2)  Autosomal dominant Noonan syndrome (BRAF, KRAS, MRAS, NRAS, PTPN11, RAF1, RIT1, SOS1, SOS2)  Autosomal dominant agenesis of corpus callosum, cardiac, ocular, and genital syndrome (CDH2)  Autosomal recessive carnitine palmitoyltransferase II deficiency (CPT2)  Autosomal dominant and autosomal recessive myofibrillar myopathy (CRYAB)  Autosomal dominant and autosomal recessive cataracts (CRYAB)  X-linked muscular dystrophy (DMD, EMD)  Autosomal recessive 3-methylglutaconic aciduria type 5 (DNAJC19)  Autosomal recessive DK1-congenital disorder of glycosylation (DOLK)  Autosomal recessive combined oxidative phosphorylation defect type 17 (ELAC2)  Autosomal recessive muscular dystrophy-dystroglycanopathy (FKRP, FKTN)  Autosomal recessive glycogen storage disease type II (GAA)  X-linked Fabry disease (GLA)  Autosomal dominant sinus node dysfunction or bradycardia (HCN4)  Autosomal dominant Costello syndrome (HRAS)  X-linked Danon disease (LAMP2)  Autosomal dominant schwannomatosis (LZTR1)  Autosomal recessive oxidative phosphorylation deficiency-10 (MTO1)  Autosomal recessive MYPN-related myopathy (MYPN)  Structural heart defects (NKX2-5, TBX20)  Autosomal recessive propionic acidemia (PCCA, PCCB)  Autosomal dominant glycogen storage-related Wolff-Parkinson-White syndrome (PRKAG2)  Autosomal dominant catecholaminergic polymorphic ventricular tachycardia (RYR2)  Autosomal recessive limb-girdle muscular dystrophy (SGCD, TCAP)</p>

Item	Description
<b>Conditions (continued)</b>	Autosomal dominant Noonan-like syndrome with loose anagen hair (SHOC2) Autosomal recessive systemic primary carnitine deficiency (SLC22A5) X-linked Barth syndrome (TFAZZIN) Autosomal recessive mitochondrial complex V deficiency (TMEM70)
<b>Interpretation</b>	All detected variants are evaluated according to American College of Medical Genetics and Genomics recommendations. Variants are classified based on known, predicted, or possible pathogenicity and reported with interpretive comments detailing their potential or known significance.
<b>Reclassification Of Variants</b>	Helix does not systematically review their variant database looking for classification changes. Helix will review the classification of previously reported variants upon request of the ordering physician/provider. Ordering physicians/providers may contact Helix Customer Support or their Dedicated Advisor and request a review of the variant classification to be performed. At the discretion of the laboratory director, the frequency of reclassification requests may be limited to once per year, no earlier than 12 months after initial variant interpretation has been performed.
<b>Variant Evaluation</b>	Variant classification is performed using the guidelines set forth by the American College of Medical Genetics and Genomics and the Association for Molecular Pathology, with modifications as suggested by domain specific Expert Panels of the Clinical genome Resource (ClinGen) when available. Variant pathogenicity is categorized as benign, likely benign, variant of uncertain significance (VUS), likely pathogenic, or pathogenic.
<b>Turnaround Time</b>	7 to 24 days
<b>Available In NY State</b>	No
<b>Test Classification</b>	This test was developed, and its performance characteristics determined, by Helix, Inc. in a manner consistent with CLIA requirements. This test has not been cleared or approved by the US Food and Drug Administration.
<b>Performing Laboratory Information</b>	CLIA Laboratory Number: 05D2117342 Laboratory Hours of Operation: Monday-Saturday (7AM-10:30PM PST) Address: 10170 Sorrento Valley Road, Suite 100, San Diego, CA 92121 Helix Customer Service: (844) 211-2070 Email: support@helix.com
<b>Regulatory Information</b>	CLIA Complexity: High Test Classification: Non-Waived/ Laboratory Developed Test
<b>CLIA Category</b>	Chemistry / Routine Chemistry

# Methods & Limitations for Helix Comprehensive Cardiomyopathy Panel



Extracted DNA is enriched for targeted regions and then sequenced using the Helix Exome+ (R) assay on an Illumina DNA sequencing system. Data is then aligned to a modified version of GRCh38 and all genes are analyzed using the MANE transcript and MANE Plus Clinical transcript, when available. Small variant calling is completed using a customized version of Sentieon's DNaseq software, augmented by a proprietary small variant caller for difficult variants. Copy number variants (CNVs) are then called using a proprietary bioinformatics pipeline based on depth analysis with a comparison to similarly sequenced samples. Interpretation is based upon guidelines published by the American College of Medical Genetics and Genomics (ACMG) and the Association for Molecular Pathology (AMP) or their modification by ClinGen Variant Curation Expert Panels when available. Interpretation is limited to the transcripts indicated on the report, +/- 10 bp into intronic regions, except as noted below. Helix variant classifications include pathogenic, likely pathogenic, variant of uncertain significance (VUS), likely benign, and benign. Variants classified as pathogenic, likely pathogenic, or VUS are included in the report. All reported variants (except for VUSs with limited evidence of pathogenicity) are confirmed through secondary manual inspection of DNA sequence data or orthogonal testing. Benign and likely benign variants are not reported but are available upon request. Risk estimations and management guidelines included in this report are based on analysis of primary literature and recommendations of applicable professional societies, and should be regarded as approximations.

Based on validation studies, this assay delivers > 99% sensitivity and specificity for single nucleotide variants and insertions and deletions (indels) up to 20 bp. Larger indels and complex variants are also reported but sensitivity may be reduced. Based on validation studies, this assay delivers > 99% sensitivity to multi-exon CNVs and > 90% sensitivity to single-exon CNVs. This test may not detect variants in challenging regions (such as short tandem repeats, homopolymer runs, and segment duplications), sub-exonic CNVs, chromosomal aneuploidy, or variants in the presence of mosaicism. Phasing will be attempted and reported, when possible. Structural rearrangements such as inversions, translocations, and gene conversions are not tested in this assay unless explicitly indicated. Additionally, deep intronic, promoter, and enhancer regions may not be covered. It is important to note that this assay cannot detect all variants known to increase disease risk, and that a negative result does not guarantee that the tested individual does not carry a rare, undetectable variant in genes analyzed. Any potential incidental findings outside of these genes and conditions will not be identified, nor reported. The results of a genetic test may be influenced by various factors, including bone marrow transplantation, blood transfusions, or in rare cases, hematolymphoid neoplasms.

## Gene Specific Notes:

*AGL*: Evaluation of chr1: 99916398 (c.4260-12A>G) will be performed. *ALPK3*: Sensitivity in exon 1 may be reduced. *BRAF*: Sensitivity in exon 1 may be reduced. *CDH2*: Sensitivity in exon 1 may be reduced. *DMD*: Evaluation of chrX:33174335 (c.31+36947G>A), chrX:31261663 (c.9225-647A>G), and chrX:31261301 (c.9225-285A>G) will be performed. *FKTN*: Evaluation of chr9:105606576 (c.648-1243G>T) will be performed. *GAA*: Evaluation of chr17:80104542 (c.-32-13T>G), chr17:80104552 (c.-32-3C>A), chr17:80104554 (c.-32-1G>C), and chr17:80108467 (c.1076-22T>G) will be performed. *GLA*: Evaluation of chrX: 101399747 (c.640-801G>A) will be performed. *MAP2K2*: Sensitivity in exon 1 may be reduced. *MYBPC3*: Evaluation of chr11:47332275-47332299 (c.3628-41\_2628-17del25), chr11:47347065 (c.906-36G>A), chr11:47346372 (c.927-2A>G), chr11:47343281 (c.1224-19G>A), chr11:47343314 (c.1224-52G>A), chr11:47343158 (c.1227-13G>A), and chr11:47340403 (c.1927+600C>T) will be performed. *PRDM16*: Analysis for exon 1 will not be performed. *PRKAG2*: Sensitivity in exon 5 may be reduced. *SLC22A5*: Evaluation of chr5:132369824 (c.-149G>A), chr5:132378362 (c.394-16T>A), and chr5:132386973 (c.825-52G>A) will be performed. *SOS2*: Sensitivity in exon 1 may be reduced. *TTN*: Analysis for exons 172 to 197 will not be performed.

## Disclaimer:

This test was developed and validated by Helix, Inc. This test has not been cleared or approved by the United States Food and Drug Administration (FDA). The Helix laboratory is accredited by the College of American Pathologists (CAP) and certified under the Clinical Laboratory Improvement Amendments (CLIA #: 05D2117342) to perform high-complexity clinical tests. This test is used for clinical purposes. It should not be regarded as investigational or for research.

# Targeted Genes & Methodology for Helix Comprehensive Cardiomyopathy Panel



The following applies to the Helix Comprehensive Cardiomyopathy Panel. Testing is performed to evaluate for the presence of variants in coding regions and extending to +/- 10 base pairs of adjacent intronic sequence on either side of the coding exons of the genes analyzed. In addition, the analysis will cover select non-coding variants, as listed below. Next-generation sequencing is performed to test for the presence of small variants and copy number variants in the genes analyzed. Confirmation of select reportable variants may be performed by alternate methodologies based on internal laboratory criteria.

This list is current from March 2024 to the present. This document is intended to highlight additional evaluations for variants of high clinical interest as well as technical limitations. For questions regarding genes, reference transcripts, or specific regions covered, contact Helix Customer Service at (844) 211-2070.

Genomic Build: GRCh38  
Catalog Number: CCMP1

Gene	Transcript	Additional Evaluations	Technical Limitations
<i>ABCC9</i>	NM_020297.4	–	–
<i>ACAD9</i>	NM_014049.5	–	–
<i>ACADVL</i>	NM_000018.4	–	–
<i>ACTC1</i>	NM_005159.5	–	–
<i>ACTN2</i>	NM_001103.4	–	–
<i>AGL</i>	NM_000642.3	Chr1: 99916398 (c.4260-12A>G)	–
<i>ALMS1</i>	NM_001378454.1	–	–
<i>ALPK3</i>	NM_020778.5	–	Sensitivity in ALPK3 exon1 may be reduced
<i>BAG3</i>	NM_004281.4	–	–
<i>BMP10</i>	NM_014482.3	–	–
<i>BRAF</i>	NM_004333.6; NM_001374258.1	–	Sensitivity to BRAF exon1 may be reduced
<i>CDH2</i>	NM_001792.5	–	Sensitivity in CDH2 exon1 may be reduced
<i>CPT2</i>	NM_000098.3	–	–
<i>CRYAB</i>	NM_001289808.2	–	–
<i>CSRP3</i>	NM_003476.5	–	–
<i>DES</i>	NM_001927.4	–	–
<i>DMD</i>	NM_004006.3	ChrX:33174335 (c.31+36947G>A) ChrX:31261663 (c.9225-647A>G) ChrX:31261301 (c.9225-285A>G)	–
<i>DNAJC19</i>	NM_145261.4	–	–
<i>DOLK</i>	NM_014908.4	–	–

# Targeted Genes & Methodology for Helix Comprehensive Cardiomyopathy Panel



Gene	Transcript	Additional Evaluations	Technical Limitations
<i>DSC2</i>	NM_024422.6	–	–
<i>DSG2</i>	NM_001943.5	–	–
<i>DSP</i>	NM_004415.4	–	–
<i>DTNA</i>	NM_001386795.1	–	–
<i>ELAC2</i>	NM_018127.7	–	–
<i>EMD</i>	NM_000117.3	–	–
<i>FHL1</i>	NM_001159699.2; NM_001159702.3	–	–
<i>FKRP</i>	NM_024301.5	–	–
<i>FKTN</i>	NM_001079802.2	Chr9:105606576 (c.648-1243G>T)	–
<i>FLNC</i>	NM_001458.5	–	–
<i>GAA</i>	NM_000152.5	Chr17:80104542 (c.-32-13T>G) Chr17:80104552 (c.-32-3C>A) Chr17:80104554 (c.-32-1G>C) Chr17:80108467 (c.1076-22T>G)	–
<i>GLA</i>	NM_000169.3	ChrX: 101399747 (c.640-801G>A)	–
<i>HCN4</i>	NM_005477.3	–	–
<i>HRAS</i>	NM_005343.4; NM_176795.5	–	–
<i>JPH2</i>	NM_020433.5	–	–
<i>JUP</i>	NM_002230.4	–	–
<i>KRAS</i>	NM_004985.5; NM_033360.4	–	–
<i>LAMP2</i>	NM_002294.3	–	–
<i>LMNA</i>	NM_170707.4; NM_005572.4	–	–
<i>LZTR1</i>	NM_006767.4	–	–
<i>MAP2K1</i>	NM_002755.4	–	–
<i>MAP2K2</i>	NM_030662.4	–	Sensitivity in MAP2K2 exon 1 may be reduced
<i>MRAS</i>	NM_001085049.3	–	–
<i>MTO1</i>	NM_012123.4	–	–

# Targeted Genes & Methodology for Helix Comprehensive Cardiomyopathy Panel



Gene	Transcript	Additional Evaluations	Technical Limitations
<i>MYBPC3</i>	NM_000256.3	Chr11:47332275-47332299 (c.3628-41_2628-17del25) Chr11:47347065 (c.906-36G>A) Chr11:47346372 (c.927-2A>G) Chr11:47343281 (c.1224-19G>A) Chr11:47343314 (c.1224-52G>A) Chr11:47343158 (c.1227-13G>A) Chr11:47340403 (c.1927+600C>T)	–
<i>MYH7</i>	NM_000257.4	–	–
<i>MYL2</i>	NM_000432.4	–	–
<i>MYL3</i>	NM_000258.3	–	–
<i>MYLK3</i>	NM_182493.3	–	–
<i>MYPN</i>	NM_032578.4	–	–
<i>NEXN</i>	NM_144573.4	–	–
<i>NKX2-5</i>	NM_004387.4	–	–
<i>NRAS</i>	NM_002524.5	–	–
<i>PCCA</i>	NM_000282.4	–	–
<i>PCCB</i>	NM_000532.5	–	–
<i>PKP2</i>	NM_001005242.3	–	–
<i>PLN</i>	NM_002667.5	–	–
<i>PPA2</i>	NM_176869.3	–	–
<i>PPCS</i>	NM_024664.4	–	–
<i>PRDM16</i>	NM_022114.4	–	Analysis for exon 1 will not be performed
<i>PRKAG2</i>	NM_016203.4	–	Sensitivity in PRKAG2 exon 5 may be reduced
<i>PTPN11</i>	NM_002834.5	–	–
<i>RAF1</i>	NM_002880.4	–	–
<i>RBM20</i>	NM_001134363.3	–	–
<i>RIT1</i>	NM_006912.6	–	–
<i>RYR2</i>	NM_001035.3	–	–
<i>SCN5A</i>	NM_000335.5; NM_001099404.2	–	–
<i>SGCD</i>	NM_000337.6	–	–
<i>SHOC2</i>	NM_007373.4	–	–

# Targeted Genes & Methodology for Helix Comprehensive Cardiomyopathy Panel



Gene	Transcript	Additional Evaluations	Technical Limitations
<i>SLC22A5</i>	NM_003060.4	Chr5:132369824 (c.-149G>A) Chr5:132378362 (c.394-16T>A) Chr5:132386973 (c.825-52G>A)	–
<i>SOS1</i>	NM_005633.4	–	–
<i>SOS2</i>	NM_006939.4	–	Sensitivity in <i>SOS2</i> exon 1 may be reduced
<i>SYNE2</i>	NM_182914.3	–	–
<i>TAFAZZIN</i>	NM_000116.5	–	–
<i>TBX20</i>	NM_001077653.2	–	–
<i>TCAP</i>	NM_003673.4	–	–
<i>TMEM43</i>	NM_024334.3	–	–
<i>TMEM70</i>	NM_017866.6	–	–
<i>TNNC1</i>	NM_003280.3	–	–
<i>TNNI3</i>	NM_000363.5	–	–
<i>TNNI3K</i>	NM_015978.3	–	–
<i>TNNT2</i>	NM_001276345.2	–	–
<i>TPM1</i>	NM_001018005.2	–	–
<i>TRIM63</i>	NM_032588.4	–	–
<i>TTN</i>	NM_001267550.2; NM_133379.5	–	Analysis for exons 172 to 197 will not be performed
<i>TTR</i>	NM_000371.4	–	–
<i>VCL</i>	NM_014000.3	–	–