

Item	Description
<b>Test Name</b>	Helix Family Variant Testing
<b>Test Type</b>	Target Analysis
<b>Catalog Number</b>	FAVT1
<b>Procedure Code</b>	H00724-6 (Helix)
<b>Test Description</b>	Helix Family Variant Testing is a targeted test to identify the presence or absence of one or more specific variants previously identified as being present in a family member. The entire gene will be evaluated and therefore additional variants determined to be pathogenic or likely pathogenic within the gene ordered will also be included in the report. Variants of uncertain significance (VUS) will not be included except in cases where the variant specified in the order is determined to be a VUS. A separate order is required for variants in separate genes.
<b>Genes Tested</b>	Based on order.
<b>Genetics Information</b>	This test utilizes next-generation sequencing to detect single nucleotide variants, insertions and deletions up to 20 bp, and copy number variants. This test includes the targeted variant or variants ordered for a given gene, along with any other pathogenic or likely pathogenic variants detected in the gene.
<b>Indications For Testing</b>	A family member of a proband who has received a genetic test result with a pathogenic or likely pathogenic variant.
<b>Clinical Descriptions</b>	Useful for diagnostic testing when a variant associated with a specific condition has been previously identified in a family member.
<b>Conditions</b>	Dependent upon the specific gene(s) ordered.
<b>Interpretation</b>	All detected variants are evaluated according to American College of Medical Genetics and Genomics recommendations. Variants are classified based on known, predicted, or possible pathogenicity; however, this test only reports pathogenic and likely pathogenic variants along with interpretive comments detailing the evidence applied towards classification. Variants of uncertain significance are not reported.
<b>Reclassification of Variants</b>	Helix does not systematically review their variant database looking for classification changes. Helix will review the classification of previously reported variants upon request of the ordering physician/provider. Ordering physicians/providers may contact Helix Customer Support or their Dedicated Advisor and request a review of the variant classification to be performed. At the discretion of the laboratory director, the frequency of reclassification requests may be limited to once per year, no earlier than 12 months after initial variant interpretation has been performed.
<b>Variant Evaluation</b>	Variant classification is performed using the guidelines set forth by the American College of Medical Genetics and Genomics and the Association for Molecular Pathology, with modifications as suggested by domain specific Expert Panels of the Clinical genome Resource (ClinGen) when available. Variant pathogenicity is categorized as benign, likely benign, variant of uncertain significance (VUS), likely pathogenic, or pathogenic.
<b>Turnaround Time</b>	7 to 24 days
<b>Available in NY State</b>	No

Item	Description
<b>Test Classification</b>	This test was developed, and its performance characteristics determined, by Helix, Inc. in a manner consistent with CLIA requirements. This test has not been cleared or approved by the US Food and Drug Administration.
<b>Performing Laboratory Information</b>	CLIA Laboratory Number: 05D2117342 Laboratory Hours of Operation: Monday-Saturday (7AM-10:30PM PST) Address: 10170 Sorrento Valley Road, Suite 100, San Diego, CA 92121 Helix Customer Service: (844) 211-2070 Email: <a href="mailto:support@helix.com">support@helix.com</a>
<b>Regulatory Information</b>	CLIA Complexity: High Test Classification: Non-Waived/ Laboratory Developed Test
<b>CLIA Category</b>	Chemistry / Routine Chemistry

Extracted DNA is enriched for targeted regions and then sequenced using the Helix Exome+ (R) assay on an Illumina DNA sequencing system. Data is then aligned to a modified version of GRCh38 and all genes are analyzed using the MANE transcript and MANE Plus Clinical transcript, when available. Small variant calling is completed using a customized version of Sentieon's DNaseq software, augmented by a proprietary small variant caller for difficult variants. Copy number variants (CNVs) are then called using a proprietary bioinformatics pipeline based on depth analysis with a comparison to similarly sequenced samples. Reportable variants in *PMS2* exons 12-15 are confirmed by PacBio long reads. The *MSH2* Boland inversion (exons 1-7) is detected by identifying discordant read-pairs spanning the presumed breakpoint. Interpretation is based upon guidelines published by the American College of Medical Genetics and Genomics (ACMG) and the Association for Molecular Pathology (AMP) or their modification by ClinGen Variant Curation Expert Panels when available. Interpretation is limited to the transcripts indicated on the report, +/- 10 bp into intronic regions, except as noted below. Helix variant classifications include pathogenic, likely pathogenic, variant of uncertain significance (VUS), likely benign, and benign. Variants classified as pathogenic, likely pathogenic, or VUS are included in the report. All reported variants (except for VUSs with limited evidence of pathogenicity) are confirmed through secondary manual inspection of DNA sequence data or orthogonal testing. Benign and likely benign variants are not reported but are available upon request. Risk estimations and management guidelines included in this report are based on analysis of primary literature and recommendations of applicable professional societies, and should be regarded as approximations.

Based on validation studies, this assay delivers > 99% sensitivity and specificity for single nucleotide variants and insertions and deletions (indels) up to 20 bp. Larger indels and complex variants are also reported but sensitivity may be reduced. Based on validation studies, this assay delivers > 99% sensitivity to multi-exon CNVs and > 90% sensitivity to single-exon CNVs. This test may not detect variants in challenging regions (such as short tandem repeats, homopolymer runs, and segment duplications), sub-exonic CNVs, chromosomal aneuploidy, or variants in the presence of mosaicism. Phasing will be attempted and reported, when possible. Structural rearrangements such as inversions, translocations, and gene conversions are not tested in this assay unless explicitly indicated. Additionally, deep intronic, promoter, and enhancer regions may not be covered. It is important to note that this assay cannot detect all variants known to increase disease risk, and that a negative result does not guarantee that the tested individual does not carry a rare, undetectable variant in genes analyzed. Any potential incidental findings outside of these genes and conditions will not be identified, nor reported. The results of a genetic test may be influenced by various factors, including bone marrow transplantation, blood transfusions, or in rare cases, hematolymphoid neoplasms.

Results do not rule out the presence of other variants related to this condition. Results may include incidental findings within the gene that are determined to be pathogenic or likely pathogenic.

#### Gene Specific Notes:

N/A

#### Disclaimer:

This test was developed and validated by Helix, Inc. This test has not been cleared or approved by the United States Food and Drug Administration (FDA). The Helix laboratory is accredited by the College of American Pathologists (CAP) and certified under the Clinical Laboratory Improvement Amendments (CLIA #: 05D2117342) to perform high-complexity clinical tests. This test is used for clinical purposes. It should not be regarded as investigational or for research.

# Targeted Genes & Methodology for Helix Family Variant Testing



The following applies to the Helix Family Variant Testing. Testing is performed to evaluate for the presence of variants in coding regions and extending to +/- 10 base pairs of adjacent intronic sequences on either side of the coding exons of the genes analyzed. In addition, the analysis will cover select non-coding variants, as listed below. Next-generation sequencing is performed to test for the presence of small variants and copy number variants in the genes analyzed. Confirmation of select reportable variants may be performed by alternate methodologies based on internal laboratory criteria.

This list is current from November 2024 to the present. This document is intended to highlight additional evaluations for variants of high clinical interest as well as technical limitations. For questions regarding genes, reference transcripts, or specific regions covered, contact Helix Customer Service at (844) 211-2070.

Genomic Build: GRCh38  
Catalog Number: FAVT1

Gene	Transcript	Additional Evaluations	Technical Limitations
<i>ABCC9</i>	NM_020297.4	–	–
<i>ACAD9</i>	NM_014049.5	–	–
<i>ACADVL</i>	NM_000018.4	–	–
<i>ACTA2</i>	NM_001613.4	–	–
<i>ACTC1</i>	NM_005159.5	–	–
<i>ACTN2</i>	NM_001103.4	–	–
<i>ADAMTS10</i>	NM_030957.4	–	–
<i>AGL</i>	NM_000642.3	Chr1: 99916398 (c.4260-12A>G)	–
<i>AIP</i>	NM_003977.4	–	–
<i>ALK</i>	NM_004304.5	–	–
<i>ALMS1</i>	NM_001378454.1	–	–
<i>ALPK3</i>	NM_020778.5	–	Sensitivity in ALPK3 exon1 may be reduced
<i>APC</i>	NM_000038.6	Includes CNV detection of Promoters 1A and 1B and sequencing of Promoter 1B	–
<i>APOB</i>	NM_000384.3	–	–
<i>ATM</i>	NM_000051.4	–	–
<i>AXIN2</i>	NM_004655.4	–	–
<i>BAG3</i>	NM_004281.4	–	–
<i>BAP1</i>	NM_004656.4	–	–
<i>BARD1</i>	NM_000465.4	–	–
<i>BGN</i>	NM_001711.6	–	–
<i>BLM</i>	NM_000057.4	–	–
<i>BMP10</i>	NM_014482.3	–	–

# Targeted Genes & Methodology for Helix Family Variant Testing



Gene	Transcript	Additional Evaluations	Technical Limitations
<i>BMPR1A</i>	NM_004329.3	Includes CNV detection in the promoter	–
<i>BRAF</i>	NM_004333.6; NM_001374258.1	–	Sensitivity to BRAF exon1 may be reduced
<i>BRCA1</i>	NM_007294.4	Sequencing analysis extends to CDS +/-20 bp	–
<i>BRCA2</i>	NM_000059.4	Sequencing analysis extends to CDS +/-20 bp	–
<i>BRIP1</i>	NM_032043.3	–	–
<i>CBS</i>	NM_000071.3	–	–
<i>CDC73</i>	NM_024529.5	–	–
<i>CDH1</i>	NM_004360.5	–	–
<i>CDH2</i>	NM_001792.5	–	Sensitivity in CDH2 exon1 may be reduced
<i>CDK4</i>	NM_000075.4	–	–
<i>CDKN1B</i>	NM_004064.5	–	–
<i>CDKN2A</i>	NM_000077.5; NM_058195.4	Includes analysis of both the p16 (p16INK4a) and p14 (p14ARF) transcripts	–
<i>CHEK2</i>	NM_007194.4	–	–
<i>COL3A1</i>	NM_000090.4	–	–
<i>COL5A1</i>	NM_000093.5	–	–
<i>COL5A2</i>	NM_000393.5	–	–
<i>CPT2</i>	NM_000098.3	–	–
<i>CRYAB</i>	NM_001289808.2	–	–
<i>CSRP3</i>	NM_003476.5	–	–
<i>CTNNA1</i>	NM_001903.5	–	–
<i>DES</i>	NM_001927.4	–	–
<i>DICER1</i>	NM_177438.3	–	–
<i>DMD</i>	NM_004006.3	ChrX:33174335 (c.31+36947G>A) ChrX:31261663 (c.9225-647A>G) ChrX:31261301 (c.9225-285A>G)	–
<i>DNAJC19</i>	NM_145261.4	–	–
<i>DOLK</i>	NM_014908.4	–	–
<i>DSC2</i>	NM_024422.6	–	–
<i>DSG2</i>	NM_001943.5	–	–

# Targeted Genes & Methodology for Helix Family Variant Testing



Gene	Transcript	Additional Evaluations	Technical Limitations
<i>DSP</i>	NM_004415.4	–	–
<i>DTNA</i>	NM_001386795.1	–	–
<i>EFEMP2</i>	NM_016938.5	–	–
<i>EGFR</i>	NM_005228.5	–	Results limited to NM_005228(EGFR): c.2369C>T (p.Thr790Met)
<i>ELAC2</i>	NM_018127.7	–	–
<i>EMD</i>	NM_000117.3	–	–
<i>EPCAM</i>	NM_002354.3	–	Results limited to CNV and limited to exons 8 and 9
<i>FBN1</i>	NM_000138.5	–	–
<i>FBN2</i>	NM_001999.4	–	–
<i>FH</i>	NM_000143.4	–	–
<i>FHL1</i>	NM_001159699.2; NM_001159702.3	–	–
<i>FKRP</i>	NM_024301.5	–	–
<i>FKTN</i>	NM_001079802.2	Chr9:105606576 (c.648-1243G>T)	–
<i>FLCN</i>	NM_144997.7	–	–
<i>FLNA</i>	NM_001110556.2	–	–
<i>FLNC</i>	NM_001458.5	–	–
<i>FOXE3</i>	NM_012186.3	–	Analysis begins at chr1:47416567 (GRCh38) and excludes the first quarter of exon 1
<i>GAA</i>	NM_000152.5	Chr17:80104542 (c.-32-13T>G) Chr17:80104552 (c.-32-3C>A) Chr17:80104554 (c.-32-1G>C) Chr17:80108467 (c.1076-22T>G)	–
<i>GLA</i>	NM_000169.3	ChrX: 101399747 (c.640-801G>A)	–
<i>GREM1</i>	NM_013372.7	–	Results limited to CNV of promoter region
<i>HCN4</i>	NM_005477.3	–	–
<i>HOXB13</i>	NM_006361.6	–	Results limited to NM_006361.6(HOXB13): c.251G>A (p.Gly84Glu)
<i>HRAS</i>	NM_005343.4; NM_176795.5	–	–
<i>JPH2</i>	NM_020433.5	–	–

# Targeted Genes & Methodology for Helix Family Variant Testing



Gene	Transcript	Additional Evaluations	Technical Limitations
<i>JUP</i>	NM_002230.4	–	–
<i>KIT</i>	NM_000222.3	–	–
<i>KRAS</i>	NM_004985.5; NM_033360.4	–	–
<i>LAMP2</i>	NM_002294.3	–	–
<i>LDLR</i>	NM_000527.5	Chr19: 11117009 (c.1845+11C>G) Chr19: 11089400 (c.-149C>A) Chr19: 11089414 (c.-135C>G) Chr19: 11089413 (c.-136C>T) Chr19: 11110640 (c.941-12G>A)	–
<i>LDLRAP1</i>	NM_015627.3	Chr1: 25564565 (c748-608G>A)	–
<i>LMNA</i>	NM_170707.4; NM_005572.4	–	–
<i>LOX</i>	NM_002317.7	–	–
<i>LZTR1</i>	NM_006767.4	–	–
<i>LZTR1</i>	NM_006767.4	–	–
<i>MAP2K1</i>	NM_002755.4	–	–
<i>MAP2K2</i>	NM_030662.4	–	Sensitivity in MAP2K2 exon 1 may be reduced
<i>MAX</i>	NM_002382.5	–	–
<i>MBD4</i>	NM_001276270.2	–	–
<i>MED12</i>	NM_005120.3	–	–
<i>MEN1</i>	NM_001370259.2	–	–
<i>MET</i>	NM_000245.4	–	–
<i>MFAP5</i>	NM_003480.4	–	–
<i>MITF</i>	NM_000248.4	–	Results limited to NM_000248.4(MITF): c.952G>A (p.Glu318Lys)
<i>MLH1</i>	NM_000249.4	Includes CNV detection in the promoter	–
<i>MRAS</i>	NM_001085049.3	–	–
<i>MSH2</i>	NM_000251.3	Includes detection of the exon 1-7 rearrangement known as the Boland Inversion and of c.942+3A>T	–
<i>MSH3</i>	NM_002439.5	–	Excludes known repeat region in MSH3 exon 1
<i>MSH6</i>	NM_000179.3	–	–
<i>MTO1</i>	NM_012123.4	–	–

# Targeted Genes & Methodology for Helix Family Variant Testing



Gene	Transcript	Additional Evaluations	Technical Limitations
<i>MUTYH</i>	NM_001048174.2; NM_001128425.2	–	–
<i>MYBPC3</i>	NM_000256.3	Chr11:47332275-47332299 (c.3628-41_2628-17del25) Chr11:47347065 (c.906-36G>A) Chr11:47346372 (c.927-2A>G) Chr11:47343281 (c.1224-19G>A) Chr11:47343314 (c.1224-52G>A) Chr11:47343158 (c.1227-13G>A) Chr11:47340403 (c.1927+600C>T)	–
<i>MYH11</i>	NM_002474.3; NM_001040113.2	–	–
<i>MYH7</i>	NM_000257.4	–	–
<i>MYL2</i>	NM_000432.4	–	–
<i>MYL3</i>	NM_000258.3	–	–
<i>MYLK</i>	NM_053025.4	–	–
<i>MYLK3</i>	NM_182493.3	–	–
<i>MYPN</i>	NM_032578.4	–	–
<i>NEXN</i>	NM_144573.4	–	–
<i>NF1</i>	NM_001042492.3	–	–
<i>NF2</i>	NM_000268.4	–	–
<i>NKX2-5</i>	NM_004387.4	–	–
<i>NOTCH1</i>	NM_017617.5	–	–
<i>NRAS</i>	NM_002524.5	–	–
<i>NTHL1</i>	NM_002528.7	–	–
<i>PALB2</i>	NM_024675.4	–	–
<i>PCCA</i>	NM_000282.4	–	–
<i>PCCB</i>	NM_000532.5	–	–
<i>PCSK9</i>	NM_174936.4	–	–
<i>PDGFRA</i>	NM_006206.6	–	–
<i>PKP2</i>	NM_001005242.3	–	–
<i>PLN</i>	NM_002667.5	–	–
<i>PLOD1</i>	NM_000302.4	–	–
<i>PMS2</i>	NM_000535.7	–	–
<i>POLD1</i>	NM_002691.4	–	CNVs not reported, and sequencing isolated to the exonuclease domain



# Targeted Genes & Methodology for Helix Family Variant Testing



Gene	Transcript	Additional Evaluations	Technical Limitations
<i>POLE</i>	NM_006231.4	–	CNVs not reported, and sequencing isolated to the exonuclease domain
<i>POT1</i>	NM_015450.3	–	–
<i>PPA2</i>	NM_176869.3	–	–
<i>PPCS</i>	NM_024664.4	–	–
<i>PRDM16</i>	NM_022114.4	–	Analysis for exon 1 will not be performed
<i>PRKAG2</i>	NM_016203.4	–	Sensitivity in PRKAG2 exon 5 may be reduced
<i>PRKAR1A</i>	NM_002734.5	–	–
<i>PRKG1</i>	NM_006258.4	–	–
<i>PTCH1</i>	NM_000264.5; NM_001083603.3	–	Sensitivity in PTCH1 exon 1 may be reduced
<i>PTEN</i>	NM_000314.8	Includes CNV detection in the promoter	–
<i>PTPN11</i>	NM_002834.5	–	–
<i>RAD51C</i>	NM_058216.3	–	–
<i>RAD51D</i>	NM_002878.4	–	–
<i>RAF1</i>	NM_002880.4	–	–
<i>RB1</i>	NM_000321.3	–	–
<i>RBM20</i>	NM_001134363.3	–	–
<i>RET</i>	NM_020975.6	–	–
<i>RIT1</i>	NM_006912.6	–	–
<i>RYR2</i>	NM_001035.3	–	–
<i>SCN5A</i>	NM_000335.5; NM_001099404.2	–	–
<i>SDHA</i>	NM_004168.4	–	CNVs not reported
<i>SDHAF2</i>	NM_017841.4	–	–
<i>SDHB</i>	NM_003000.3	–	–
<i>SDHC</i>	NM_003001.5	–	–
<i>SDHD</i>	NM_003002.4	–	–
<i>SGCD</i>	NM_000337.6	–	–
<i>SHOC2</i>	NM_007373.4	–	–
<i>SKI</i>	NM_003036.4	–	–
<i>SLC22A5</i>	NM_003060.4	Chr5:132369824 (c.-149G>A) Chr5:132378362 (c.394-16T>A) Chr5:132386973 (c.825-52G>A)	–

# Targeted Genes & Methodology for Helix Family Variant Testing



Gene	Transcript	Additional Evaluations	Technical Limitations
<i>SLC2A10</i>	NM_030777.4	–	–
<i>SMAD2</i>	NM_005901.6	–	–
<i>SMAD3</i>	NM_005902.4	–	–
<i>SMAD4</i>	NM_005359.6	–	–
<i>SMAD4</i>	NM_005359.6	–	–
<i>SMARCA4</i>	NM_003072.5; NM_001387283.1	–	–
<i>SMARCB1</i>	NM_003073.5	–	–
<i>SMARCE1</i>	NM_003079.5	–	–
<i>SOS1</i>	NM_005633.4	–	–
<i>SOS2</i>	NM_006939.4	–	Sensitivity in <i>SOS2</i> exon 1 may be reduced
<i>STK11</i>	NM_000455.5	–	Sensitivity in <i>STK11</i> exon 3 may be reduced
<i>SUFU</i>	NM_016169.4	–	–
<i>SYNE2</i>	NM_182914.3	–	–
<i>TAFAZZIN</i>	NM_000116.5	–	–
<i>TBX20</i>	NM_001077653.2	–	–
<i>TCAP</i>	NM_003673.4	–	–
<i>TGFB2</i>	NM_003238.6	–	–
<i>TGFB3</i>	NM_003239.5	–	–
<i>TGFBR1</i>	NM_004612.4	–	Analysis for exon 1 will not be performed
<i>TGFBR2</i>	NM_003242.6	–	–
<i>TMEM127</i>	NM_017849.4	–	–
<i>TMEM43</i>	NM_024334.3	–	–
<i>TMEM70</i>	NM_017866.6	–	–
<i>TNNC1</i>	NM_003280.3	–	–
<i>TNNI3</i>	NM_000363.5	–	–
<i>TNNI3K</i>	NM_015978.3	–	–
<i>TNNT2</i>	NM_001276345.2	–	–
<i>TP53</i>	NM_000546.6	Includes CNV detection in the promoter	–
<i>TPM1</i>	NM_001018005.2	–	–
<i>TRIM63</i>	NM_032588.4	–	–

# Targeted Genes & Methodology for Helix Family Variant Testing



Gene	Transcript	Additional Evaluations	Technical Limitations
<i>TSC1</i>	NM_000368.5	–	Sensitivity in TSC1 exon 21 may be reduced
<i>TSC2</i>	NM_000548.5	–	–
<i>TTN</i>	NM_001267550.2; NM_133379.5	–	Analysis for exons 172 to 197 will not be performed
<i>TTR</i>	NM_000371.4	–	–
<i>VCL</i>	NM_014000.3	–	–
<i>VHL</i>	NM_000551.4	–	Excludes coverage of cryptic exon E1' (chr3:10142758-10143009)