Helix STAT Breast Cancer Panel



Item	Description		
Test Name	Helix STAT Breast Cancer Panel		
Test Type	Hereditary Cancer		
Catalog Number	SBRC1		
Procedure Code	H00825-3 (Helix)		
Test Description	This STAT panel evaluates 13 genes that have an established primary association with hereditary predisposition to breast cancer and delivers expedited results to enable accelerated clinical intervention.		
Genes Tested	ATM, BARD1, BRCA1, BRCA2, CDH1, CHEK2, NF1, PALB2, PTEN, RAD51C, RAD51D, STK11, TP53		
Genetics Information	This test utilizes next-generation sequencing to detect single nucleotide variants, insertions and deletions up to 20 bp, and copy number variants in genes associated with hereditary forms of breast cancer.		
Indications For Testing	A personal and/or family history suggestive of a hereditary form of breast cancer predisposition.		
Clinical Descriptions	Hereditary predisposition to breast cancer refers to the increased likelihood of developing adult-onset breast cancer due to genetic factors that are passed down in families. The most well-known genes are <i>BRCA1</i> and <i>BRCA2</i> , although there are others, including the genes on this panel. Individuals with a pathogenic variant may also have an increased risk of other cancers such as ovarian and prostate, depending on the affected gene. The genes on this panel were specifically selected for their established association with breast cancer and the availability of published medical management and surveillance recommendations. These recommendations may include increased cancer screening and preventive surgery for early-detection and prevention. Identification of a pathogenic variant also helps identify at-risk family members, who can pursue genetic testing and preventive measures. The genes on this panel are associated with conditions that have autosomal dominant and/or autosomal recessive inheritance. Note that some of these genes may also be associated with other unrelated conditions; this means that when undergoing this test, there is a possibility of incidentally detecting carrier status for, or predisposition to, one of these conditions.		
Conditions	Hereditary breast and ovarian cancer syndrome (BRCA1 and BRCA2) Hereditary diffuse gastric cancer (CDH1) PALB2-related cancer susceptibility (PALB2) PTEN hamartoma tumor syndrome (PTEN) Peutz-Jeghers syndrome (STK11) Li-Fraumeni syndrome (TP53) ATM-related cancer susceptibility (ATM) BARD1-related cancer susceptibility (BARD1) CHEK2-related cancer susceptibility (CHEK2) Neurofibromatosis (NF1) RAD51C-related cancer susceptibility (RAD51C) RAD51D-related cancer susceptibility (RAD51D)		

Helix STAT Breast Cancer Panel



Item	Description		
	In addition, one or more of the genes on this panel are associated with other conditions for what a predisposition to, or carrier status of, may be identified: Autosomal recessive Fanconi anemia (BRCA1, BRCA2, PALB2 and RAD51C) Autosomal recessive ataxia-telangiectasia (ATM)		
Interpretation	All detected variants are evaluated according to American College of Medical Genetics and Genomics recommendations. Variants are classified based on known, predicted, or possible pathogenicity and reported with interpretive comments detailing their potential or known significance.		
Reclassification Of Variants	Helix reviews variant classifications annually when they arise in routine processes and upon request from providers. The timing of re-review depends on clinical risk. Providers can request a variant re-review by contacting Helix Customer Support. If a classification by Helix is updated, Helix identifies affected past patients and issues revised reports. Updated results are communicated to providers prior to results being uploaded to the EHR, and patients are notified through the EHR patient portal.		
Variant Evaluation	Variant classification is performed using the guidelines set forth by the American College of Medical Genetics and Genomics and the Association for Molecular Pathology, with modifications as suggested by domain specific Expert Panels of the Clinical genome Resource (ClinGen) when available. Variant pathogenicity is categorized as benign, likely benign, variant of uncertain significance (VUS), likely pathogenic, or pathogenic.		
Turnaround Time - Standard	Typically 7 to 11 days		
Turnaround Time - Requery (SOQO®)	Typically ≤ 2 days		
Available In NY State	Yes		
Test Classification	This test was developed, and its performance characteristics determined, by Helix, Inc. in a manner consistent with CLIA requirements. This test has not been cleared or approved by the US Food and Drug Administration.		
Performing Laboratory Information	CLIA Laboratory Number: 05D2117342 Laboratory Hours of Operation: Monday-Saturday (7AM-10:30PM PST) Address: 10170 Sorrento Valley Road, Suite 100, San Diego, CA 92121 Helix Customer Service: (844) 211-2070 Email: support@helix.com		
Regulatory Information	Test Classification: Non-Waived/ Laboratory Developed Test		
CLIA Category	Chemistry / Routine Chemistry		

Methods & Limitations for Helix STAT Breast Cancer Panel



Extracted DNA is enriched for targeted regions and then sequenced using the Helix Exome+ (R) assay on an Illumina DNA sequencing system. Data is then aligned to a modified version of GRCh38 and all genes are analyzed using the MANE transcript and MANE Plus Clinical transcript, when available. Small variant calling is completed using a customized version of Sentieon's DNAseq software, augmented by a proprietary small variant caller for difficult variants. Copy number variants (CNVs) are then called using a proprietary bioinformatics pipeline based on depth analysis with a comparison to similarly sequenced samples. Interpretation is based upon guidelines published by the American College of Medical Genetics and Genomics (ACMG) and the Association for Molecular Pathology (AMP) or their modification by ClinGen Variant Curation Expert Panels when available. Interpretation is limited to the transcripts indicated on the report, +/- 10 bp into intronic regions, except as noted below. Helix variant classifications include pathogenic, likely pathogenic, variant of uncertain significance (VUS), likely benign, and benign. Variants classified as pathogenic, likely pathogenic, or VUS are included in the report. All reported variants (except for VUSs with limited evidence of pathogenicity) are confirmed through secondary manual inspection of DNA sequence data or orthogonal testing. Benign and likely benign variants are not reported but are available upon request. Risk estimations and management guidelines included in this report are based on analysis of primary literature and recommendations of applicable professional societies, and should be regarded as approximations.

Based on validation studies, this assay delivers > 99% sensitivity and specificity for single nucleotide variants and insertions and deletions (indels) up to 20 bp. Larger indels and complex variants are also reported but sensitivity may be reduced. Based on validation studies, this assay delivers > 99% sensitivity to multi-exon CNVs and > 90% sensitivity to single-exon CNVs. This test may not detect variants in challenging regions (such as short tandem repeats, homopolymer runs, and segment duplications), sub-exonic CNVs, chromosomal aneuploidy, or variants in the presence of mosaicism. Phasing will be attempted and reported, when possible. Structural rearrangements such as inversions, translocations, and gene conversions are not tested in this assay unless explicitly indicated. Additionally, deep intronic, promoter, and enhancer regions may not be covered. It is important to note that this assay cannot detect all variants known to increase disease risk, and that a negative result does not guarantee that the tested individual does not carry a rare, undetectable variant in genes analyzed. Any potential incidental findings outside of these genes and conditions will not be identified, nor reported. The results of a genetic test may be influenced by various factors, including bone marrow transplantation, blood transfusions, or in rare cases, hematolymphoid neoplasms.

Gene Specific Notes:

BRCA1: sequencing analysis extends to CDS +/-20 bp; *BRCA2*: sequencing analysis extends to CDS +/-20 bp; *PTEN*: analysis includes CNV of the promoter; *STK11*: sensitivity of exon 3 analysis may be reduced; *TP53*: analysis includes CNV of the promoter.

Disclaimer:

This test was developed and validated by Helix, Inc. This test has not been cleared or approved by the United States Food and Drug Administration (FDA). The Helix laboratory is accredited by the College of American Pathologists (CAP) and certified under the Clinical Laboratory Improvement Amendments (CLIA #: 05D2117342) to perform high-complexity clinical tests. This test is used for clinical purposes. It should not be regarded as investigational or for research.

Targeted Genes & Methodology for Helix STAT Breast Cancer Panel



The following applies to the Helix STAT Breast Cancer Panel. Testing is performed to evaluate for the presence of variants in coding regions and extending to +/- 10 base pairs of adjacent intronic sequences on either side of the coding exons of the genes analyzed. In addition, the analysis will cover select non-coding variants, as listed below. Next-generation sequencing is performed to test for the presence of small variants and copy number variants in the genes analyzed. Confirmation of select reportable variants may be performed by alternate methodologies based on internal laboratory criteria.

This list is current from July 2025 to the present. This document is intended to highlight additional evaluations for variants of high clinical interest as well as technical limitations. For questions regarding genes, reference transcripts, or specific regions covered, contact Helix Customer Service at (844) 211-2070.

Genomic Build: GRCh38 Catalog Number: SBRC1

Gene	Transcript	Additional Evaluations	Technical Limitations
ATM	NM_000051.4	-	-
BARD1	NM_000465.4	-	-
BRCA1	NM_007294.4	Sequencing analysis extends to CDS +/-20 bp	_
BRCA2	NM_000059.4	Sequencing analysis extends to CDS +/-20 bp	-
CDH1	NM_004360.5	-	-
CHEK2	NM_007194.4	-	-
NF1	NM_001042492.3	-	-
PALB2	NM_024675.4	-	-
PTEN	NM_000314.8	Includes CNV detection in the promoter	-
RAD51C	NM_058216.3	-	-
RAD51D	NM_002878.4	-	-
STK11	NM_000455.5	_	Sensitivity in STK11 exon 3 may be reduced
TP53	NM_000546.6	Includes CNV detection in the promoter	_