



## Overview

### Test Description

Helix Family Variant Testing is a targeted test to identify the presence or absence of one or more specific variants previously identified as being present in a family member. The entire gene will be evaluated for each gene included in the order, based on the analytical and reportable range described below, and therefore additional variants determined to be pathogenic or likely pathogenic within the genes ordered will also be included in the report. Variants of uncertain significance (VUS) will not be included except in cases where the variant specified in the order is determined to be a VUS.

### Genetics Information

This test utilizes next-generation sequencing to detect single nucleotide variants, insertions and deletions up to 20 bp, and copy number variants. This test includes the targeted variant or variants ordered for a given gene included in the order, along with any other pathogenic or likely pathogenic variants detected in the gene(s).

### Indications For Testing

A family member of a proband who has received a genetic test result with a test result with an identified variant within the analytical and reportable range of this test.

### Clinical Descriptions

Useful for diagnostic testing when the patient is at increased risk to carry a variant associated with a specific condition that has been previously identified in a family member.

### Conditions

Dependent upon the specific gene(s) ordered.

## Test Details

### Genes Tested

*Based on order.*

### Methods & Limitations

Extracted DNA is enriched for targeted regions and then sequenced using the Helix Exome+ (R) assay on an Illumina DNA sequencing system. Data is then aligned to a modified version of GRCh38 and all genes are analyzed using the MANE transcript and MANE Plus Clinical transcript, when available. Small variant calling is completed using a customized version of Sentieon's DNaseq software, augmented by a proprietary small variant caller for difficult variants. Copy number variants (CNVs) are then called using a proprietary bioinformatics pipeline based on depth analysis with a comparison to similarly sequenced samples. Reportable variants in *PMS2* exons 11-15 are confirmed by PacBio long reads. Both the *MSH2* Boland inversion (exons 1-7) and the *BRCA2* Alu insertion are detected by identifying discordant read-pairs spanning the presumed breakpoint. Interpretation is based upon guidelines published by the American College of Medical Genetics and Genomics (ACMG) and the Association for Molecular Pathology (AMP) or their modification by ClinGen Variant Curation Expert Panels when available. Interpretation is limited to the transcripts indicated on the report, +/- 10 bp into intronic regions, except as noted below. Helix variant classifications include pathogenic, likely pathogenic, variant of uncertain significance (VUS), likely benign, and benign. Variants classified as pathogenic, likely pathogenic, or VUS if the requested variant is classified as such, are included in the report. All reported variants (except for VUSs with limited evidence of pathogenicity) are confirmed through secondary manual inspection of DNA sequence data or orthogonal testing. Benign and likely benign variants are not reported but are available upon request. Risk estimations and management guidelines included in this report are based on analysis of primary literature and recommendations of applicable professional societies, and should be regarded as approximations.

Based on validation studies, this assay delivers > 99% sensitivity and specificity for single nucleotide variants and insertions and deletions (indels) up to 20 bp. Larger indels and complex variants are also reported but sensitivity may be reduced. Based on validation studies, this assay delivers > 99% sensitivity to multi-exon CNVs and > 90% sensitivity to single-exon CNVs. This test may not detect variants in challenging regions (such as short tandem repeats, homopolymer runs, and segmental duplications), sub-exonic CNVs, chromosomal aneuploidy, or variants in the presence of mosaicism. Phasing will be attempted and reported, when possible. Structural rearrangements such as inversions, translocations, complex rearrangements, and gene conversions are not tested in this assay unless explicitly indicated. Additionally, deep intronic, promoter, and enhancer regions may not be covered. It is important to note that this assay cannot detect all variants known to increase disease risk, and that a negative result does not guarantee that the tested individual does not carry a rare, undetectable variant in the genes analyzed. Any potential incidental findings outside of these genes and conditions will not be identified, nor reported. The results of a genetic test may be influenced by various factors, including bone marrow transplantation, blood transfusions, or in rare cases, hematolymphoid neoplasms.

Results do not rule out the presence of other variants related to this condition. Results may include incidental findings within the gene that are determined to be pathogenic or likely pathogenic.



**Gene Specific Notes:**

Gene specific notes vary by ordered gene. For details, visit [helix.com/clinical-genomics/targeted-analyses/family-variant-testing](https://helix.com/clinical-genomics/targeted-analyses/family-variant-testing). For questions, contact [clinicalsupport@helix.com](mailto:clinicalsupport@helix.com).

**Disclaimer**

This test was developed and validated by Helix, Inc. This test has not been cleared or approved by the United States Food and Drug Administration (FDA). The Helix laboratory is accredited by the College of American Pathologists (CAP) and certified under the Clinical Laboratory Improvement Amendments (CLIA #: 05D2117342) to perform high-complexity clinical tests. This test is used for clinical purposes. It should not be regarded as investigational use only or for research use only.

**Interpretation**

All detected variants are evaluated according to American College of Medical Genetics and Genomics recommendations. Variants are classified based on known, predicted, or possible pathogenicity; however, this test only reports pathogenic and likely pathogenic variants along with interpretive comments detailing the evidence applied towards classification. Variants of uncertain significance are not reported unless specifically requested.

**Reclassification**

Helix reviews variant classifications annually when they arise in routine processes and upon request from providers. The timing of re-review depends on clinical risk. Providers can request a variant re-review by contacting Helix Customer Support. If a classification by Helix is updated, Helix identifies affected past patients and issues revised reports. Updated results are communicated to providers prior to results being uploaded to the EHR, and patients are notified through the EHR patient portal.

**Variant Evaluation**

Variant classification is performed using the guidelines set forth by the American College of Medical Genetics and Genomics and the Association for Molecular Pathology, with modifications as suggested by domain specific Expert Panels of the Clinical genome Resource (ClinGen) when available. Variant pathogenicity is categorized as benign, likely benign, variant of uncertain significance (VUS), likely pathogenic, or pathogenic.

**Laboratory Details**

**Turnaround Time - Standard**

Typically 6 to 21 days

**Turnaround Time - Requery (SOQO<sup>®</sup>)**

Typically ≤ 5 days

**Available in NY State**

Yes

**Procedure Code**

H00824-6 (Helix)



Laboratory Details

Test Classification

This test was developed, and its performance characteristics determined, by Helix, Inc. in a manner consistent with CLIA requirements. This test has not been cleared or approved by the US Food and Drug Administration.

Regulatory Information

- CLIA Complexity: High
- Test Classification: Non-Waived/ Laboratory Developed Test

CLIA Category

Chemistry / Routine Chemistry

Performing Laboratory Information

- CLIA Laboratory Number: 05D2117342
- Laboratory Hours of Operation: Monday-Saturday (7AM-10:30PM PST)
- Address: 10170 Sorrento Valley Road, Suite 100, San Diego, CA 92121
- Helix Customer Service: (844) 211-2070
- Email: [support@helix.com](mailto:support@helix.com)



## Helix Family Variant Testing

FAVT2 | Diagnostic Test

The following applies to the Helix Family Variant Testing. Testing is performed to evaluate for the presence of variants in coding regions and extending to +/- 10 base pairs of adjacent intronic sequences on either side of the coding exons of the genes analyzed. In addition, the analysis will cover select non-coding variants, as listed below. Next-generation sequencing is performed to test for the presence of small variants and copy number variants in the genes analyzed. Confirmation of select reportable variants may be performed by alternate methodologies based on internal laboratory criteria.

This list is current from April 2026 to the present. This document is intended to highlight additional evaluations for variants of high clinical interest as well as technical limitations. For questions regarding genes, reference transcripts, or specific regions covered, contact Helix Customer Service at (844) 211-2070.

Genomic Build: GRCh38

Gene	Transcript	Special Gene Notes
<i>ABCA1</i>	NM_005502.4	Analysis for exon 40 will not be performed
<i>ABCC9</i>	NM_020297.4	–
<i>ABCD1</i>	NM_000033.4	–
<i>ABCG5</i>	NM_022436.3	–
<i>ABCG8</i>	NM_022437.3	–
<i>ACAD9</i>	NM_014049.5	–
<i>ACADVL</i>	NM_000018.4	–
<i>ACTA2</i>	NM_001613.4	–
<i>ACTC1</i>	NM_005159.5	–
<i>ACTN2</i>	NM_001103.4	–
<i>ACVRL1</i>	NM_000020.3	–
<i>ADAMTS10</i>	NM_030957.4	–
<i>AGL</i>	NM_000642.3	Chr1: 99916398 (c.4260-12A>G)
<i>AIP</i>	NM_003977.4	–
<i>ALK</i>	NM_004304.5	–
<i>ALMS1</i>	NM_001378454.1	–
<i>ALPK3</i>	NM_020778.5	Sensitivity in <i>ALPK3</i> exon1 may be reduced
<i>ANGPTL3</i>	NM_014495.4	–
<i>ANK2</i>	NM_001148.6	–
<i>APC</i>	NM_000038.6	Includes CNV detection of Promoters 1A and 1B and sequencing of Promoter 1B
<i>APOA1</i>	NM_000039.3	–



# Helix Family Variant Testing

FAVT2 | Diagnostic Test

Gene	Transcript	Special Gene Notes
<i>APOA5</i>	NM_001371904.1	–
<i>APOB</i>	NM_000384.3	–
<i>APOC2</i>	NM_000483.5	–
<i>ATM</i>	NM_000051.4	–
<i>ATP7B</i>	NM_000053.4	–
<i>AXIN2</i>	NM_004655.4	–
<i>BAG3</i>	NM_004281.4	–
<i>BAP1</i>	NM_004656.4	–
<i>BARD1</i>	NM_000465.4	–
<i>BGN</i>	NM_001711.6	–
<i>BLM</i>	NM_000057.4	–
<i>BMP10</i>	NM_014482.3	–
<i>BMPR1A</i>	NM_004329.3	Includes CNV detection in the promoter
<i>BRAF</i>	NM_004333.6; NM_001374258.1	Sensitivity to <i>BRAF</i> exon1 may be reduced
<i>BRCA1</i>	NM_007294.4	Sequencing analysis extends to CDS +/-20 bp
<i>BRCA2</i>	NM_000059.4	Includes c.156_157insAlu and Sequencing analysis extends to CDS +/-20 bp
<i>BRIP1</i>	NM_032043.3	–
<i>BTD</i>	NM_001370658.1	–
<i>CACNA1C</i>	NM_000719.7; NM_001167623.2	–
<i>CACNA1D</i>	NM_001128840.3; NM_000720.4	–
<i>CACNA1S</i>	NM_000069.3	Sensitivity to exon 91 may be reduced
<i>CALM1</i>	NM_006888.6	–
<i>CALM2</i>	NM_001743.6	–
<i>CALM3</i>	NM_005184.4	–
<i>CASQ2</i>	NM_001232.4	–
<i>CAV3</i>	NM_033337.3	–
<i>CBS</i>	NM_000071.3	–



# Helix Family Variant Testing

FAVT2 | Diagnostic Test

Gene	Transcript	Special Gene Notes
<i>CDC73</i>	NM_024529.5	–
<i>CDH1</i>	NM_004360.5	–
<i>CDH2</i>	NM_001792.5	Sensitivity in <i>CDH2</i> exon1 may be reduced
<i>CDK4</i>	NM_000075.4	–
<i>CDKN1B</i>	NM_004064.5	–
<i>CDKN2A</i>	NM_000077.5; NM_058195.4	Includes analysis of both the p16 (p16INK4a) and p14 (p14ARF) transcripts
<i>CEBPA</i>	NM_004364.5	–
<i>CETP</i>	NM_000078.3	–
<i>CHEK2</i>	NM_007194.4	–
<i>COL3A1</i>	NM_000090.4	–
<i>COL5A1</i>	NM_000093.5	–
<i>COL5A2</i>	NM_000393.5	–
<i>CPT2</i>	NM_000098.3	–
<i>CREB3L3</i>	NM_032607.3	–
<i>CRYAB</i>	NM_001289808.2	–
<i>CSRP3</i>	NM_003476.5	–
<i>CTNNA1</i>	NM_001903.5	–
<i>CYP27A1</i>	NM_000784.4	–
<i>DDX41</i>	NM_016222.4	–
<i>DES</i>	NM_001927.4	–
<i>DICER1</i>	NM_177438.3	–
<i>DMD</i>	NM_004006.3	ChrX:33174335 (c.31+36947G>A) ChrX:31261663 (c.9225-647A>G) ChrX:31261301 (c.9225-285A>G)
<i>DNAJC19</i>	NM_145261.4	–
<i>DOLK</i>	NM_014908.4	–
<i>DSC2</i>	NM_024422.6	–
<i>DSG2</i>	NM_001943.5	–
<i>DSP</i>	NM_004415.4	–
<i>DTNA</i>	NM_001386795.1	–



Gene	Transcript	Special Gene Notes
<i>EFEMP2</i>	NM_016938.5	–
<i>EGFR</i>	NM_005228.5	Results limited to NM_005228( <i>EGFR</i> ): c.2369C>T (p.Thr790Met)
<i>ELAC2</i>	NM_018127.7	–
<i>EMD</i>	NM_000117.3	–
<i>ENG</i>	NM_001114753.3	–
<i>EPCAM</i>	NM_002354.3	Results limited to CNV and limited to exons 8 and 9
<i>ETV6</i>	NM_001987.5	–
<i>F2</i>	NM_000506.5	Results limited to c.*97G>A
<i>F5</i>	NM_000130.5	Results limited to c.1601G>A (p.Arg534Gln)
<i>FBN1</i>	NM_000138.5	–
<i>FBN2</i>	NM_001999.4	–
<i>FH</i>	NM_000143.4	–
<i>FHL1</i>	NM_001159699.2; NM_001159702.3	–
<i>FKRP</i>	NM_024301.5	–
<i>FKTN</i>	NM_001079802.2	Chr9:105606576 (c.648-1243G>T)
<i>FLCN</i>	NM_144997.7	–
<i>FLNA</i>	NM_001110556.2	–
<i>FLNC</i>	NM_001458.5	–
<i>FOXE3</i>	NM_012186.3	Analysis begins at chr1:47416567 (GRCh38) and excludes the first quarter of exon 1
<i>GAA</i>	NM_000152.5	Chr17:80104542 (c.-32-13T>G) Chr17:80104552 (c.-32-3C>A) Chr17:80104554 (c.-32-1G>C) Chr17:80108467 (c.1076-22T>G)
<i>GATA2</i>	NM_032638.5; NM_001145661.2	–
<i>GLA</i>	NM_000169.3	ChrX: 101399747 (c.640-801G>A)
<i>GPD1</i>	NM_005276.4	–
<i>GPIHBP1</i>	NM_178172.6	–
<i>GREM1</i>	NM_013372.7	Results limited to CNV of promoter region
<i>HCN4</i>	NM_005477.3	–



# Helix Family Variant Testing

FAVT2 | Diagnostic Test

Gene	Transcript	Special Gene Notes
<i>HFE</i>	NM_000410.4	Results limited to c.845G>A (p.Cys282Tyr) and c.187C>G (p.His63Asp)
<i>HNF1A</i>	NM_000545.8	–
<i>HOXB13</i>	NM_006361.6	Results limited to NM_006361.6( <i>HOXB13</i> ): c.251G>A (p.Gly84Glu)
<i>HRAS</i>	NM_005343.4; NM_176795.5	–
<i>JPH2</i>	NM_020433.5	–
<i>JUP</i>	NM_002230.4	–
<i>KCNE1</i>	NM_000219.6	–
<i>KCNE2</i>	NM_172201.2	–
<i>KCNH2</i>	NM_000238.4	Chr7:150958048-150958065 (c.910_916+11del)
<i>KCNJ2</i>	NM_000891.3	–
<i>KCNQ1</i>	NM_000218.3	Sensitivity in <i>KCNQ1</i> exon 1 may be reduced Chr11:2461715 (c.386+16231G>A) Chr11:2585210-2585211 (c.1033-1_1117dup)
<i>KIT</i>	NM_000222.3	–
<i>KRAS</i>	NM_004985.5; NM_033360.4	–
<i>LAMP2</i>	NM_002294.3	–
<i>LCAT</i>	NM_000229.2	–
<i>LDLR</i>	NM_000527.5	–
<i>LDLRAP1</i>	NM_015627.3	–
<i>LIPA</i>	NM_000235.4	–
<i>LIPG</i>	NM_006033.4	–
<i>LMF1</i>	NM_022773.4	–
<i>LMNA</i>	NM_170707.4; NM_005572.4	–
<i>LOX</i>	NM_002317.7	–
<i>LPL</i>	NM_000237.3	–
<i>LRP6</i>	NM_002336.3	–
<i>LZTR1</i>	NM_006767.4	–
<i>MAP2K1</i>	NM_002755.4	–
<i>MAP2K2</i>	NM_030662.4	Sensitivity in <i>MAP2K2</i> exon 1 may be reduced



Gene	Transcript	Special Gene Notes
MAX	NM_002382.5	–
MBD4	NM_001276270.2	–
MED12	NM_005120.3	–
MEN1	NM_001370259.2	–
MET	NM_000245.4	–
MFAP5	NM_003480.4	–
MITF	NM_000248.4	Results limited to NM_000248.4(MITF): c.952G>A (p.Glu318Lys)
MLH1	NM_000249.4	Includes CNV detection in the promoter
MRAS	NM_001085049.3	–
MSH2	NM_000251.3	Includes detection of the exon 1-7 rearrangement known as the Boland Inversion and of c.942+3A>T
MSH3	NM_002439.5	Excludes known repeat region in MSH3 exon 1
MSH6	NM_000179.3	–
MTO1	NM_012123.4	–
MTTP	NM_001386140.1	–
MUTYH	NM_001048174.2; NM_001128425.2	–
MYBPC3	NM_000256.3	Chr11:47332275-47332299 (c.3628-41_2628-17del25) Chr11:47347065 (c.906-36G>A) Chr11:47346372 (c.927-2A>G) Chr11:47343281 (c.1224-19G>A) Chr11:47343314 (c.1224-52G>A) Chr11:47343158 (c.1227-13G>A) Chr11:47340403 (c.1927+600C>T)
MYH11	NM_002474.3; NM_001040113.2	–
MYH7	NM_000257.4	–
MYL2	NM_000432.4	–
MYL3	NM_000258.3	–
MYL4	NM_002476.2	–
MYLK	NM_053025.4	–
MYLK3	NM_182493.3	–
MYPN	NM_032578.4	–



# Helix Family Variant Testing

FAVT2 | Diagnostic Test

Gene	Transcript	Special Gene Notes
<i>NEXN</i>	NM_144573.4	–
<i>NF1</i>	NM_001042492.3	–
<i>NF2</i>	NM_000268.4	–
<i>NKX2-5</i>	NM_004387.4	–
<i>NOTCH1</i>	NM_017617.5	–
<i>NRAS</i>	NM_002524.5	–
<i>NTHL1</i>	NM_002528.7	–
<i>OTC</i>	NM_000531.6	–
<i>PALB2</i>	NM_024675.4	–
<i>PCCA</i>	NM_000282.4	–
<i>PCCB</i>	NM_000532.5	–
<i>PCSK9</i>	NM_174936.4	–
<i>PDGFRA</i>	NM_006206.6	–
<i>PHOX2B</i>	NM_003924.4	Excludes exon 3, which encompasses the alanine repeat region
<i>PKP2</i>	NM_001005242.3	–
<i>PLN</i>	NM_002667.5	–
<i>PLOD1</i>	NM_000302.4	–
<i>PMS2</i>	NM_000535.7	–
<i>PNPLA2</i>	NM_020376.4	–
<i>POLD1</i>	NM_002691.4	CNVs not reported, and sequencing isolated to the exonuclease domain
<i>POLE</i>	NM_006231.4	CNVs not reported, and sequencing isolated to the exonuclease domain
<i>POT1</i>	NM_015450.3	–
<i>PPA2</i>	NM_176869.3	–
<i>PPARG</i>	NM_138711.6	–
<i>PPCS</i>	NM_024664.4	–
<i>PRDM16</i>	NM_022114.4	Analysis for exon 1 will not be performed
<i>PRKAG2</i>	NM_016203.4	Sensitivity in <i>PRKAG2</i> exon 5 may be reduced
<i>PRKAR1A</i>	NM_002734.5	–
<i>PRKG1</i>	NM_006258.4	–
<i>PTCH1</i>	NM_000264.5; NM_001083603.3	Sensitivity in <i>PTCH1</i> exon 1 may be reduced



# Helix Family Variant Testing

FAVT2 | Diagnostic Test

Gene	Transcript	Special Gene Notes
<i>PTEN</i>	NM_000314.8	Includes CNV detection in the promoter
<i>PTPN11</i>	NM_002834.5	–
<i>RAD51C</i>	NM_058216.3	–
<i>RAD51D</i>	NM_002878.4	–
<i>RAF1</i>	NM_002880.4	–
<i>RB1</i>	NM_000321.3	–
<i>RBM20</i>	NM_001134363.3	–
<i>RET</i>	NM_020975.6	–
<i>RIT1</i>	NM_006912.6	–
<i>RPE65</i>	NM_000329.3	–
<i>RPS20</i>	NM_001023.4	–
<i>RUNX1</i>	NM_001754.5	–
<i>RYR1</i>	NM_000540.3	–
<i>RYR2</i>	NM_001035.3	–
<i>SAR1B</i>	NM_016103.4	Sensitivity in <i>SAR1B</i> exon 4 may be reduced
<i>SCN5A</i>	NM_000335.5; NM_001099404.2	–
<i>SDHA</i>	NM_004168.4	CNVs not reported
<i>SDHAF2</i>	NM_017841.4	–
<i>SDHB</i>	NM_003000.3	–
<i>SDHC</i>	NM_003001.5	–
<i>SDHD</i>	NM_003002.4	–
<i>SERPINA1</i>	NM_000295.5, NM_001127701.1	Results limited to c.1096G>A (Z allele), c.863A>T (S allele)
<i>SGCD</i>	NM_000337.6	–
<i>SHOC2</i>	NM_007373.4	–
<i>SKI</i>	NM_003036.4	–
<i>SLC22A5</i>	NM_003060.4	Chr5:132369824 (c.-149G>A) Chr5:132378362 (c.394-16T>A) Chr5:132386973 (c.825-52G>A)
<i>SLC2A10</i>	NM_030777.4	–



# Helix Family Variant Testing

FAVT2 | Diagnostic Test

Gene	Transcript	Special Gene Notes
<i>SLC4A3</i>	NM_005070.4	–
<i>SMAD2</i>	NM_005901.6	–
<i>SMAD3</i>	NM_005902.4	–
<i>SMAD4</i>	NM_005359.6	–
<i>SMARCA4</i>	NM_003072.5; NM_001387283.1	–
<i>SMARCB1</i>	NM_003073.5	–
<i>SMARCE1</i>	NM_003079.5	–
<i>SOS1</i>	NM_005633.4	–
<i>SOS2</i>	NM_006939.4	Sensitivity in <i>SOS2</i> exon 1 may be reduced
<i>STAC3</i>	NM_145064.3	Sensitivity to exon 7 may be reduced
<i>STK11</i>	NM_000455.5	Sensitivity in <i>STK11</i> exon 3 may be reduced
<i>SUFU</i>	NM_016169.4	–
<i>SYNE2</i>	NM_182914.3	–
<i>TAFAZZIN</i>	NM_000116.5	–
<i>TBX20</i>	NM_001077653.2	–
<i>TCAP</i>	NM_003673.4	–
<i>TECRL</i>	NM_001010874.5	–
<i>TGFB2</i>	NM_003238.6	–
<i>TGFB3</i>	NM_003239.5	–
<i>TGFBR1</i>	NM_004612.4	Analysis for exon 1 will not be performed
<i>TGFBR2</i>	NM_003242.6	–
<i>TMEM127</i>	NM_017849.4	–
<i>TMEM43</i>	NM_024334.3	–
<i>TMEM70</i>	NM_017866.6	–
<i>TNNC1</i>	NM_003280.3	–
<i>TNNI3</i>	NM_000363.5	–
<i>TNNI3K</i>	NM_015978.3	–
<i>TNNT2</i>	NM_001276345.2	–
<i>TP53</i>	NM_000546.6	Includes CNV detection in the promoter



Gene	Transcript	Special Gene Notes
<i>TPM1</i>	NM_001018005.2	–
<i>TRDN</i>	NM_006073.4	Chr6:123636725 (c.22+29A>G)
<i>TRIM63</i>	NM_032588.4	–
<i>TSC1</i>	NM_000368.5	Sensitivity in <i>TSC1</i> exon 21 may be reduced
<i>TSC2</i>	NM_000548.5	–
<i>TTN</i>	NM_001267550.2	Analysis for exons 172 to 197 will not be performed
<i>TTR</i>	NM_000371.4	–
<i>VCL</i>	NM_014000.3	–
<i>VHL</i>	NM_000551.4	Excludes coverage of cryptic exon E1' (chr3:10142758-10143009)
<i>WT1</i>	NM_024426.6	–